

For Research Use Only

OX40/TNFRSF4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84247-3-PBS



Basic Information

Catalog Number:

84247-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003327

GeneID (NCBI):

7293

UNIPROT ID:

P43489

Full Name:

tumor necrosis factor receptor
superfamily, member 4

Calculated MW:

29 kDa

Purification Method:

Protein A purification

CloneNo.:

241398G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01138-1, CD134/OX40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84247-2-PBS. Detection antibody: 84247-3-PBS. Standard: Eg1261. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01138-3, CD134/OX40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84247-1-PBS. Detection antibody: 84247-3-PBS. Standard: Eg1261. Range: 0.313-40 ng/mL.