

For Research Use Only

CD134/OX40 Recombinant antibody

Catalog Number: 84247-4-RR



Basic Information

Catalog Number:

84247-4-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_003327

GeneID (NCBI):

7293

UNIPROT ID:

P43489

Full Name:

tumor necrosis factor receptor
superfamily, member 4

Calculated MW:

29 kDa

Observed MW:

45-50 kDa

Purification Method:

Protein A purification

CloneNo.:

241398B4

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MJ cells,

Background Information

CD134 (OX40) is a member of the TNFR-superfamily of receptors. Predominantly expressed on activated T cells, CD134 is activated by its cognate ligand CD134L (OX40L) and functions as a T cell co-stimulatory molecule. CD134-CD134L interactions have been proposed as a potential therapeutic target for treating autoimmunity (PMID: 26215166). CD134 can interact with TRAF2, TRAF3, and TRAF5.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84247-4-RR (OX40/TNFRSF4 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84247-4-RR against Human CD134/OX40 were performed. The affinity constant is 5.21 nM.