

For Research Use Only

IL3RA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84256-2-PBS



Basic Information

Catalog Number:

84256-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC035407

GeneID (NCBI):

3563

ENSEMBL Gene ID:

ENSG00000185291

UNIPROT ID:

P26951

Full Name:

interleukin 3 receptor, alpha (low affinity)

Calculated MW:

378 aa, 43 kDa

Purification Method:

Protein A purification

CloneNo.:

241267C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00176-3, IL3RA/CD123 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84256-2-PBS. Detection antibody: 84256-6-PBS. Standard: RP02336. Range: 0.313-40 ng/mL.