

For Research Use Only

BTLA/CD272 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84262-1-PBS



Basic Information

Catalog Number: 84262-1-PBS	GenBank Accession Number: BC107091	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 151888	CloneNo.: 241318G6
Source: Rabbit	UNIPROT ID: Q7Z6A9	
Isotype: IgG	Full Name: B and T Lymphocyte associated	
	Calculated MW: 289 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01168-3, BTLA/CD272 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84262-3-PBS. Detection antibody: 84262-1-PBS. Standard: RP02326. Range: 0.078-10 ng/mL.