

For Research Use Only

CD22 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84265-3-PBS



Basic Information

Catalog Number:

84265-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC109306

GeneID (NCBI):

933

UNIPROT ID:

P20273

Full Name:

CD22 molecule

Calculated MW:

847 aa, 95 kDa

Purification Method:

Protein A purification

CloneNo.:

241314A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01173-1, CD22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84265-2-PBS. Detection antibody: 84265-3-PBS. Standard: RPO2324. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01173-2, CD22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84265-1-PBS. Detection antibody: 84265-3-PBS. Standard: RPO2324. Range: 0.313-40 ng/mL