

For Research Use Only

TOR1B Recombinant antibody, PBS Only (Capture)

Catalog Number: 84267-1-PBS



Basic Information

Catalog Number: 84267-1-PBS	GenBank Accession Number: BC015578	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 27348	CloneNo.: 241574F4
Source: Rabbit	UNIPROT ID: O14657	
Isotype: IgG	Full Name: torsin family 1, member B (torsin B)	
Immunogen Catalog Number: AG16641	Calculated MW: 336 aa, 38 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01175-1, TOR1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84267-1-PBS. Detection antibody: 84267-3-PBS. Standard: Ag16641. Range: 0.313-40 ng/mL.