

For Research Use Only

IVNS1ABP Recombinant antibody, PBS Only (Capture)

Catalog Number: 84268-2-PBS



Basic Information

Catalog Number: 84268-2-PBS	GenBank Accession Number: BC067739	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10625	CloneNo.: 241587H4
Source: Rabbit	UNIPROT ID: Q9Y6Y0	
Isotype: IgG	Full Name: influenza virus NS1A binding protein	
Immunogen Catalog Number: AG6412	Calculated MW: 72 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01182-2, IVNS1ABP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84268-2-PBS. Detection antibody: 84268-1-PBS. Standard: Ag6412. Range: 1.25-40 ng/mL.



Cytometric bead array standard curve of MP01182-3, IVNS1ABP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84268-2-PBS. Detection antibody: 84268-3-PBS. Standard: Ag6412. Range: 1.25-40 ng/mL.