

For Research Use Only

Rat Il1 beta Recombinant antibody, PBS Only (Detector)

Catalog Number: 84274-4-PBS



Basic Information

Catalog Number: 84274-4-PBS	GenBank Accession Number: NM_031512	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 24494	CloneNo.: 241001A1
Source: Rabbit	UNIPROT ID: Q63264	
Isotype: IgG	Full Name: interleukin 1 beta	
	Calculated MW: 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01210-2, RAT IL-1 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84274-2-PBS. Detection antibody: 84274-4-PBS. Standard: Eg31371. Range: 0.313-40 ng/mL.