

For Research Use Only

ZNF217 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84276-2-PBS



Basic Information

Catalog Number: 84276-2-PBS	GenBank Accession Number: BC113427	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7764	CloneNo.: 241217E5
Source: Rabbit	UNIPROT ID: O75362	
Isotype: IgG	Full Name: zinc finger protein 217	
Immunogen Catalog Number: AG19453	Calculated MW: 1048 aa, 115 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01357-1, Human ZNF217 Recombinant Matched Antibody Pair - PBS only. 84276-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag19453. 84276-2-PBS was HRP conjugated as the detection antibody. Range: 6.25-400 ng/mL

The mean ZNF217 concentration was determined to be 885 ng/mL in HepG2 cell extract based on a 1.3 mg/mL extract load.