

For Research Use Only

NME3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84278-1-PBS



Basic Information

Catalog Number: 84278-1-PBS	GenBank Accession Number: BC000250	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4832	CloneNo.: 241585D6
Source: Rabbit	UNIPROT ID: Q13232	
Isotype: IgG	Full Name: non-metastatic cells 3, protein expressed in	
Immunogen Catalog Number: AG7284	Calculated MW: 19 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01200-1, NME3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84278-2-PBS. Detection antibody: 84278-1-PBS. Standard: Ag7284. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01200-2, NME3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84278-5-PBS. Detection antibody: 84278-1-PBS. Standard: Ag7284. Range: 0.313-40 ng/mL