

For Research Use Only

FAT1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84279-1-PBS



Basic Information

Catalog Number: 84279-1-PBS	GenBank Accession Number: GeneID (NCBI): 2195	Purification Method: Protein A purification
Size: 1 mg/ml	UNIPROT ID: Q14517	CloneNo.: 241602H6
Source: Rabbit	Full Name: FAT tumor suppressor homolog 1 (Drosophila)	
Isotype: IgG		
Immunogen Catalog Number: AG33957		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01192-1, FAT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84279-1-PBS. Detection antibody: 84279-2-PBS. Standard: Ag33957. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01192-2, FAT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84279-1-PBS. Detection antibody: 84279-4-PBS. Standard: Ag33957. Range: 0.625-80 ng/mL