

For Research Use Only

# FAT1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84279-4-PBS



## Basic Information

Catalog Number: 84279-4-PBS	GenBank Accession Number: GeneID (NCBI): 2195	Purification Method: Protein A purification
Size: 1 mg/ml	UNIPROT ID: Q14517	CloneNo.: 241602G4
Source: Rabbit	Full Name: FAT tumor suppressor homolog 1 (Drosophila)	
Isotype: IgG		
Immunogen Catalog Number: AG33957		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01192-2, FAT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84279-1-PBS. Detection antibody: 84279-4-PBS. Standard: Ag33957. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01192-3, FAT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84279-4-PBS. Detection antibody: 84279-3-PBS. Standard: Ag33957. Range: 0.625-80 ng/mL