

For Research Use Only

TRAIL Recombinant antibody, PBS Only (Detector)

Catalog Number: 84280-3-PBS



Basic Information

Catalog Number: 84280-3-PBS	GenBank Accession Number: BC032722	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 8743	CloneNo.: 241419G6
Source: Rabbit	UNIPROT ID: P50591	
Isotype: IgG	Full Name: tumor necrosis factor (ligand) superfamily, member 10	
	Calculated MW: 281 aa, 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01166-2, TRAIL Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84280-2-PBS. Detection antibody: 84280-3-PBS. Standard: RP02339. Range: 0.313-40 ng/mL