

For Research Use Only

MLPH Recombinant antibody, PBS Only (Capture)

Catalog Number: 84282-3-PBS



Basic Information

Catalog Number: 84282-3-PBS	GenBank Accession Number: BC001653	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 79083	CloneNo.: 241621H2
Source: Rabbit	UNIPROT ID: Q9BV36	
Isotype: IgG	Full Name: melanophilin	
Immunogen Catalog Number: AG32931	Calculated MW: 66 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01196-2, MLPH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84282-3-PBS. Detection antibody: 84282-5-PBS. Standard: Ag32931. Range: 0.313-40 ng/mL