

For Research Use Only

# NOG Recombinant antibody, PBS Only (Capture)

Catalog Number: 84283-1-PBS



## Basic Information

Catalog Number: 84283-1-PBS	GenBank Accession Number: BC034027	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9241	CloneNo.: 241511C9
Source: Rabbit	UNIPROT ID: Q13253	
Isotype: IgG	Full Name: noggin	
Immunogen Catalog Number: AG32755	Calculated MW: 26 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01174-1, Noggin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84283-1-PBS. Detection antibody: 84283-2-PBS. Standard: Ag32755. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01174-2, Noggin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84283-1-PBS. Detection antibody: 84283-3-PBS. Standard: Ag32755. Range: 0.313-40 ng/mL