

For Research Use Only

NOG Recombinant antibody, PBS Only (Capture)

Catalog Number: 84283-6-PBS



Basic Information

Catalog Number: 84283-6-PBS	GenBank Accession Number: BC034027	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9241	CloneNo.: 241511E5
Source: Rabbit	UNIPROT ID: Q13253	
Isotype: IgG	Full Name: noggin	
Immunogen Catalog Number: AG32755	Calculated MW: 26 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01174-4, Human Noggin Recombinant Matched Antibody Pair - PBS only. 84283-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag32755. 84283-3-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL.



The mean Noggin concentration was determined to be 79.96 ng/mL in mouse brain tissue extract based on a 2.0 mg/mL extract load.