

For Research Use Only

TMEM98 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84284-1-PBS



Basic Information

Catalog Number: 84284-1-PBS	GenBank Accession Number: BC000526	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 26022	CloneNo.: 241581B9
Source: Rabbit	UNIPROT ID: Q9Y2Y6	
Isotype: IgG	Full Name: transmembrane protein 98	
Immunogen Catalog Number: AG6492	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01195-1, TMEM98 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84284-1-PBS. Detection antibody: 84284-4-PBS. Standard: Ag6492. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01195-2, TMEM98 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84284-1-PBS. Detection antibody: 84284-3-PBS. Standard: Ag6492. Range: 0.313-40 ng/mL.



Bi-layer interferometry (BLI) kinetic assay of 84284-1-PBS against Human TMEM98 was performed. The affinity constant is below 1 pM.