

For Research Use Only

TMEM98 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84284-3-PBS



Basic Information

| | | |
|--|---|---|
| Catalog Number: 84284-3-PBS | GenBank Accession Number: BC000526 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 26022 | CloneNo.: 241581A2 |
| Source: Rabbit | UNIPROT ID: Q9Y2Y6 | |
| Isotype: IgG | Full Name: transmembrane protein 98 | |
| Immunogen Catalog Number: AG6492 | Calculated MW: 25 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01195-2, TMEM98 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84284-1-PBS. Detection antibody: 84284-3-PBS. Standard: Ag6492. Range:0.313-40 ng/mL.



Bilayer interferometry (BLI) kinetic assay of 84284-3-PBS against Human TMEM98 was performed. The affinity constant is 0.30 nM.