

For Research Use Only

# FLT3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84286-3-PBS



## Basic Information

Catalog Number:

84286-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC036028

GeneID (NCBI):

2322

UNIPROT ID:

P36888

Full Name:

fms-related tyrosine kinase 3

Calculated MW:

993 aa, 113 kDa

Purification Method:

Protein A purification

CloneNo.:

241413E8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01184-2, FLT3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84286-3-PBS. Detection antibody: 84286-1-PBS. Standard: RPO2337. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01184-3, FLT3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84286-3-PBS. Detection antibody: 84286-4-PBS. Standard: RPO2337. Range: 0.625-40 ng/mL