

For Research Use Only

ADAM17 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84292-1-PBS



Basic Information

Catalog Number: 84292-1-PBS	GenBank Accession Number: BC136783	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6868	CloneNo.: 241305D1
Source: Rabbit	UNIPROT ID: P78536	
Isotype: IgG	Full Name: ADAM metallopeptidase domain 17	
Immunogen Catalog Number: AG32418	Calculated MW: 824 aa, 93 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01169-1, ADAM17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84292-2-PBS. Detection antibody: 84292-1-PBS. Standard: Ag32418. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP01169-2, ADAM17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84292-3-PBS. Detection antibody: 84292-1-PBS. Standard: Ag32418. Range: 0.625-40 ng/mL.