

For Research Use Only

GSDMC Recombinant antibody, PBS Only (Capture)

Catalog Number: 84294-2-PBS



Basic Information

Catalog Number: 84294-2-PBS	GenBank Accession Number: NM_031415	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 56169	CloneNo.: 241489F11
Source: Rabbit	UNIPROT ID: Q9BYG8	
Isotype: IgG	Full Name: gasdermin C	
Immunogen Catalog Number: AG33319	Calculated MW: 58kd	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01172-1, GSDMC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84294-2-PBS. Detection antibody: 84294-1-PBS. Standard: Ag33319. Range: 0.625-80 ng/mL