

For Research Use Only

# P glycoprotein Recombinant antibody, PBS Only (Capture)

Catalog Number:84295-1-PBS



## Basic Information

Catalog Number:

84295-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC130424

GeneID (NCBI):

5243

UNIPROT ID:

P08183

Full Name:

ATP-binding cassette, sub-family B  
(MDR/TAP), member 1

Calculated MW:

1280 aa, 142 kDa

Purification Method:

Protein A purification

CloneNo.:

241512E5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01176-1, P glycoprotein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84295-1-PBS. Detection antibody: 84295-2-PBS. Standard: Eg2366. Range: 0.313-40 ng/mL.