

For Research Use Only

P glycoprotein Recombinant antibody, PBS Only (Detector)

Catalog Number: 84295-2-PBS



Basic Information

Catalog Number: 84295-2-PBS	GenBank Accession Number: BC130424	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5243	CloneNo.: 241512B1
Source: Rabbit	UNIPROT ID: P08183	
Isotype: IgG	Full Name: ATP-binding cassette, sub-family B (MDR/TAP), member 1	
	Calculated MW: 1280 aa, 142 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01176-1, P glycoprotein Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84295-1-PBS. Detection antibody: 84295-2-PBS. Standard: Eg2366. Range: 0.313-40 ng/mL.