

For Research Use Only

Mouse Atox1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84305-2-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 84305-2-PBS | GenBank Accession Number: NM_009720 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 11927 | CloneNo.: 241566C2 |
| Source: Rabbit | UNIPROT ID: O08997 | |
| Isotype: IgG | Full Name: ATX1 (antioxidant protein 1) homolog 1 (yeast) | |
| Immunogen Catalog Number: AG35290 | Calculated MW: 7 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01193-1, MOUSE Atox1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84305-1-PBS. Detection antibody: 84305-2-PBS. Standard: Ag35290. Range: 0.313-40 ng/mL.