

For Research Use Only

FTSJ3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84311-4-PBS



Basic Information

Catalog Number: 84311-4-PBS	GenBank Accession Number: BC036710	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 117246	CloneNo.: 241507F6
Source: Rabbit	UNIPROT ID: Q8IY81	
Isotype: IgG	Full Name: FtsJ homolog 3 (E. coli)	
Immunogen Catalog Number: AG9835	Calculated MW: 847 aa, 97 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01213-2, FTSJ3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84311-4-PBS. Detection antibody: 84311-2-PBS. Standard: Ag9835. Range: 0.156-20 ng/mL