

For Research Use Only

MCTS1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84312-4-PBS



Basic Information

Catalog Number: 84312-4-PBS	GenBank Accession Number: BC001013	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 28985	CloneNo.: 241576H2
Source: Rabbit	UNIPROT ID: Q9ULC4	
Isotype: IgG	Full Name: malignant T cell amplified sequence	
Immunogen Catalog Number: AG6968	Calculated MW: 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only


For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01214-1, MCTS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84312-3-PBS. Detection antibody: 84312-4-PBS. Standard: Ag6968. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01214-3, MCTS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84312-1-PBS. Detection antibody: 84312-4-PBS. Standard: Ag6968. Range: 0.313-40 ng/mL