

For Research Use Only

Mouse 4-1BB/TNFRSF9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84315-1-PBS



Basic Information

Catalog Number: 84315-1-PBS	GenBank Accession Number: NM_001077509.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 21942	CloneNo.: 241647D8
Source: Rabbit	UNIPROT ID: P20334	
Isotype: IgG	Full Name: tumor necrosis factor receptor superfamily, member 9	
	Calculated MW: 28kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01216-2, MOUSE 4-1BB/TNFRSF9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84315-1-PBS. Detection antibody: 84315-3-PBS. Standard: Eg1627. Range: 0.078-10 ng/mL