

For Research Use Only

C14orf133 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:84320-1-PBS



Basic Information

Catalog Number: 84320-1-PBS	GenBank Accession Number: BC015054	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 63894	CloneNo.: 241616A6
Source: Rabbit	UNIPROT ID: Q9H9C1	
Isotype: IgG	Full Name: chromosome 14 open reading frame 133	
Immunogen Catalog Number: AG14762	Calculated MW: 493 aa, 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01194-1, VIPAS39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84320-1-PBS. Detection antibody: 84320-3-PBS. Standard: Ag14762. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01194-3, VIPAS39 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84320-2-PBS. Detection antibody: 84320-1-PBS. Standard: Ag14762. Range: 0.781-100 ng/mL