

For Research Use Only

ICAM3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84335-4-PBS



Basic Information

Catalog Number: 84335-4-PBS	GenBank Accession Number: BC058903	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3385	CloneNo.: 241418G8
Source: Rabbit	UNIPROT ID: P32942	
Isotype: IgG	Full Name: intercellular adhesion molecule 3	
	Calculated MW: 547 aa, 60 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01223-3, Human ICAM3/CD50 Recombinant Matched Antibody Pair - PBS only. 84335-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard RPO2342. 84335-2-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL



Serum of twelve individual healthy human donors was measured. The human ICAM3/CD50 concentration of detected samples was determined to be 100.57 ng/mL with a range of 80.46 - 116.80 ng/mL.