

For Research Use Only

ZNF148 Recombinant antibody, PBS Only (Capture)

Catalog Number:84338-4-PBS



Basic Information

Catalog Number:

84338-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35742

GenBank Accession Number:

BC035591

GeneID (NCBI):

7707

UNIPROT ID:

Q9UQR1

Full Name:

zinc finger protein 148

Purification Method:

Protein A purification

CloneNo.:

241632F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01253-3, ZNF148 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84338-4-PBS. Detection antibody: 84338-5-PBS. Standard: Ag35742. Range: 0.313-40 ng/mL