

For Research Use Only

ZNF148 Recombinant antibody, PBS Only

Catalog Number: 84338-6-PBS



Basic Information

Catalog Number:

84338-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35742

GenBank Accession Number:

BC035591

GeneID (NCBI):

7707

UNIPROT ID:

Q9UQR1

Full Name:

zinc finger protein 148

Purification Method:

Protein A purification

CloneNo.:

241632A7

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using ZNF148 antibody (84338-6-RR, Clone: 241632A7) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 84338-6-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 84338-6-RR against Human ZNF148 were performed. The affinity constant is 2.71 nM.