

For Research Use Only

# Mouse CD24 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84342-3-PBS



## Basic Information

Catalog Number:

84342-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_009846.2

GeneID (NCBI):

12484

UNIPROT ID:

P24807

Full Name:

CD24a antigen

Calculated MW:

8kDa

Purification Method:

Protein A purification

CloneNo.:

241691F5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01237-1, MOUSE CD24 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84342-1-PBS. Detection antibody: 84342-3-PBS. Standard: Eg2104. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01237-3, MOUSE CD24 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84342-4-PBS. Detection antibody: 84342-3-PBS. Standard: Eg2104. Range: 0.781-100 ng/mL