

For Research Use Only

Alpha-2-Macroglobulin Recombinant antibody, PBS Only (Capture)

Catalog Number: 84347-1-PBS



Basic Information

Catalog Number: 84347-1-PBS	GenBank Accession Number: NM_000014.6	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2	CloneNo.: 241649D11
Source: Rabbit	UNIPROT ID: P01023	
Isotype: IgG	Full Name: alpha-2-macroglobulin	
	Calculated MW: 163kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01224-1, Alpha-2-macroglobulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84347-1-PBS. Detection antibody: 84347-2-PBS. Standard: null. Range: null.



Cytometric bead array standard curve of MP01224-3, Alpha-2-macroglobulin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84347-1-PBS. Detection antibody: 84347-4-PBS. Standard: Eg1198. Range: 0.469-60 ng/mL.