

For Research Use Only

Mouse CD2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84360-3-PBS



Basic Information

Catalog Number: 84360-3-PBS	GenBank Accession Number: NM_013486.2	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 12481	CloneNo.: 241548D5
Source: Rabbit	UNIPROT ID: P08920	
Isotype: IgG	Full Name: CD2 antigen	
	Calculated MW: 38kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01249-3, MOUSE CD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84360-5-PBS. Detection antibody: 84360-3-PBS. Standard: Eg1403. Range: 0.313-40 ng/mL.