

For Research Use Only

Placental lactogen Recombinant antibody, PBS Only (Capture)

Catalog Number: 84361-4-PBS



Basic Information

Catalog Number: 84361-4-PBS	GenBank Accession Number: BC002717	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1442	CloneNo.: 241625F12
Source: Rabbit	UNIPROT ID: PODML3	
Isotype: IgG	Full Name: chorionic somatomammotropin hormone 1 (placental lactogen)	
Immunogen Catalog Number: AG0201	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01252-2, CSH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84361-4-PBS. Detection antibody: 84361-2-PBS. Standard: Ag0201. Range: 0.781-100 ng/mL.