

For Research Use Only

# IGFBP6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84364-4-PBS



## Basic Information

Catalog Number:

84364-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC011708

GeneID (NCBI):

3489

UNIPROT ID:

P24592

Full Name:

insulin-like growth factor binding protein 6

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

241724G5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01229-2, IGFBP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84364-3-PBS. Detection antibody: 84364-4-PBS. Standard: Eg0087. Range: 0.313-40 ng/mL