

For Research Use Only

THAP11 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84373-2-PBS



Basic Information

Catalog Number: 84373-2-PBS	GenBank Accession Number: BC012182	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57215	CloneNo.: 241590B1
Source: Rabbit	UNIPROT ID: Q96EK4	
Isotype: IgG	Full Name: THAP domain containing 11	
Immunogen Catalog Number: AG19342	Calculated MW: 314 aa, 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP01244-1, THAP11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84373-5-PBS. Detection antibody: 84373-2-PBS. Standard: Ag19342. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01244-3, THAP11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84373-1-PBS. Detection antibody: 84373-2-PBS. Standard: null. Range: null.



Cytometric bead array standard curve of MP01244-3, THAP11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84373-1-PBS. Detection antibody: 84373-2-PBS. Standard: Ag19342. Range: 0.781-100 ng/mL