

For Research Use Only

HAPLN1 Recombinant antibody, PBS Only (Capture)

Catalog Number:84374-5-PBS



Basic Information

Catalog Number: 84374-5-PBS	GenBank Accession Number: BC057808	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1404	CloneNo.: 241588H5
Source: Rabbit	UNIPROT ID: P10915	
Isotype: IgG	Full Name: hyaluronan and proteoglycan link protein 1	
Immunogen Catalog Number: AG25221	Calculated MW: 40 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01245-3, Human HAPLN1 Recombinant Matched Antibody Pair - PBS only. 84374-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag25221. 84374-2-PBS was HRP conjugated as the detection antibody. Range: 3.13-200 ng/mL.



The mean HAPLN1 concentration was determined to be 1,727.92 ng/mL in human placenta cell extract based on a 1.70 mg/mL extract load.