

For Research Use Only

# CYP2R1 Recombinant antibody

Catalog Number: 84378-2-RR



## Basic Information

<b>Catalog Number:</b> 84378-2-RR	<b>GenBank Accession Number:</b> BC104909	<b>Purification Method:</b> Protein A purification
<b>Size:</b>	<b>GeneID (NCBI):</b> 120227	<b>CloneNo.:</b> 241072E2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q6VVX0	
<b>Isotype:</b> IgG	<b>Full Name:</b> cytochrome P450, family 2, subfamily R, polypeptide 1	
<b>Immunogen Catalog Number:</b> AG12873	<b>Calculated MW:</b> 501 aa, 57 kDa	

## Applications

**Tested Applications:**  
ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data