

For Research Use Only

CYP2R1 Recombinant antibody

Catalog Number: 84378-4-RR



Basic Information

Catalog Number: 84378-4-RR	GenBank Accession Number: BC104909	Purification Method: Protein A purification
Size:	GeneID (NCBI): 120227	CloneNo.: 241072F5
Source: Rabbit	UNIPROT ID: Q6VVX0	
Isotype: IgG	Full Name: cytochrome P450, family 2, subfamily R, polypeptide 1	
Immunogen Catalog Number: AG12873	Calculated MW: 501 aa, 57 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data