

For Research Use Only

PDGFRA Recombinant antibody, PBS Only (Detector)

Catalog Number: 84383-1-PBS



Basic Information

Catalog Number: 84383-1-PBS	GenBank Accession Number: BC015186	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5156	CloneNo.: 241658F12
Source: Rabbit	UNIPROT ID: P16234	
Isotype: IgG	Full Name: platelet-derived growth factor receptor, alpha polypeptide	
	Calculated MW: 123 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01256-3, PDGFRA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84383-4-PBS. Detection antibody: 84383-1-PBS. Standard: RP00003. Range: 0.781-100 ng/mL