

For Research Use Only

# FCAR Recombinant antibody, PBS Only (Capture)

Catalog Number: 84384-4-PBS



## Basic Information

Catalog Number:

84384-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC027953

GeneID (NCBI):

2204

UNIPROT ID:

P24071

Full Name:

Fc fragment of IgA, receptor for

Calculated MW:

287 aa, 32 kDa

Purification Method:

Protein A purification

CloneNo.:

241426D7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01257-3, CD89 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84384-4-PBS. Detection antibody: 84384-1-PBS. Standard: RP02338. Range: 0.156-10 ng/mL