

For Research Use Only

CD98/SLC3A2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84387-1-PBS



Basic Information

Catalog Number:

84387-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002394.5

GeneID (NCBI):

6520

UNIPROT ID:

P08195-1

Full Name:

solute carrier family 3 (activators of dibasic and neutral amino acid transport), member 2

Calculated MW:

68kDa

Purification Method:

Protein A purification

CloneNo.:

241755C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01259-2, CD98/SLC3A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84387-1-PBS. Detection antibody: 84387-4-PBS. Standard: Eg1948. Range: 0.156-20 ng/mL.



Cytometric bead array standard curve of MP01259-3, CD98/SLC3A2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84387-1-PBS. Detection antibody: 84387-2-PBS. Standard: Eg1948. Range: 0.156-20 ng/mL.