

For Research Use Only

Mouse PDGFR beta Recombinant antibody, PBS Only (Capture)

Catalog Number: 84388-3-PBS



Basic Information

Catalog Number:

84388-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008809.2

GeneID (NCBI):

18596

UNIPROT ID:

P05622-1

Full Name:

platelet derived growth factor receptor, beta polypeptide

Calculated MW:

123kDa

Purification Method:

Protein A purification

CloneNo.:

241754E1

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP01260-1, Mouse PDGFR beta Recombinant Matched Antibody Pair - PBS only. 84388-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1533. 84388-1-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



Serum of six mice was measured. The mouse PDGFR beta concentration of detected samples was determined to be 34.36 ng/mL with a range of 25.40 - 43.60 ng/mL.



Cytometric bead array standard curve of MP01260-1, MOUSE PDGFR beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84388-3-PBS. Detection antibody: 84388-1-PBS. Standard: Eg1533. Range: 0.625-80 ng/mL