

For Research Use Only

IL12RB1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84391-2-PBS



Basic Information

Catalog Number: 84391-2-PBS	GenBank Accession Number: BC029121	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3594	CloneNo.: 241660E3
Source: Rabbit	UNIPROT ID: P42701	
Isotype: IgG	Full Name: interleukin 12 receptor, beta 1	
	Calculated MW: 662 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01262-1, IL12RB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84391-5-PBS. Detection antibody: 84391-2-PBS. Standard: RP00005. Range: 0.781-50 ng/mL.