

For Research Use Only

SLC28A3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84393-1-PBS



Basic Information

Catalog Number:

84393-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12758

GenBank Accession Number:

BC093821

GeneID (NCBI):

64078

UNIPROT ID:

Q9HAS3

Full Name:

solute carrier family 28 (sodium-coupled nucleoside transporter), member 3

Calculated MW:

691 aa, 77 kDa

Purification Method:

Protein A purification

CloneNo.:

241615F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01264-1, SLC28A3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84393-2-PBS. Detection antibody: 84393-1-PBS. Standard: Ag12758. Range: 0.781-100 ng/mL.