

For Research Use Only

TCEB3B Recombinant antibody, PBS Only (Detector)

Catalog Number: 84400-1-PBS



Basic Information

Catalog Number:

84400-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10595

GenBank Accession Number:

BC036022

GeneID (NCBI):

51224

UNIPROT ID:

Q8IYF1

Full Name:

transcription elongation factor B
polypeptide 3B (elongin A2)

Calculated MW:

753 aa, 84 kDa

Purification Method:

Protein A purification

CloneNo.:

241497E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01271-1, TCEB3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84400-2-PBS. Detection antibody: 84400-1-PBS. Standard: Ag10595. Range: 0.781-100 ng/mL.