

For Research Use Only

NeuN Recombinant antibody, PBS Only (Detector)

Catalog Number: 84401-3-PBS



Basic Information

Catalog Number:

84401-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25689

GenBank Accession Number:

NM_001082575

GeneID (NCBI):

146713

UNIPROT ID:

A6NFN3

Full Name:

hexaribonucleotide binding protein 3

Purification Method:

Protein A purification

CloneNo.:

241461G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01273-1, NeuN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84401-2-PBS. Detection antibody: 84401-3-PBS. Standard: Ag25689. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01273-2, NeuN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84401-1-PBS. Detection antibody: 84401-3-PBS. Standard: Ag25689. Range: 0.156-20 ng/mL