

For Research Use Only

CD40/TNFRSF5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84406-3-PBS



Basic Information

Catalog Number:

84406-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001250.6

GeneID (NCBI):

958

ENSEMBL Gene ID:

ENSG00000101017

UNIPROT ID:

P25942-1

Full Name:

CD40 molecule, TNF receptor
superfamily member 5

Calculated MW:

31kDa

Purification Method:

Protein A purification

CloneNo.:

241760C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01294-2, CD40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84406-3-PBS. Detection antibody: 84406-4-PBS. Standard: Eg1797. Range: 0.625-80 ng/mL