

For Research Use Only

# CD40/TNFRSF5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84406-4-PBS



## Basic Information

Catalog Number:

84406-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001250.6

GeneID (NCBI):

958

ENSEMBL Gene ID:

ENSG00000101017

UNIPROT ID:

P25942-1

Full Name:

CD40 molecule, TNF receptor  
superfamily member 5

Calculated MW:

31kDa

Purification Method:

Protein A purification

CloneNo.:

241760C10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01294-2, CD40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84406-3-PBS. Detection antibody: 84406-4-PBS. Standard: Eg1797. Range: 0.625-80 ng/mL